

## Caring Mug



- The Caring Mug features a wide stable base with two large contoured and angled handles
- It is designed to be held in a comfortable position for drinking and is useful for people with tremor or who cannot sit upright. The mug is microwaveable
- Weight: 176 g

### Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Eating and drinking

**Category:** Low risk

Click [here](#) for more information

### NDIS

**Item Code:**

05\_150900319\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Unit Weight	Capacity
DLK276510	Clear Mug	176 g	300 ml