

## Leckey Totstander



- Totstander is designed to suit children at early stages of standing
- Often used as an assessment tool
- Available with adjustable pelvic belt, chest strap, knee cups and sandals
- Highly portable
- User Standing Height: 750 - 900 mm
- Chest Pad Height: 470 - 640 mm



**ndis**

**Item Code:**

05\_053603010\_0103\_1\_2

Click [here](#) for more information

## Warranty Information

- 2 Years

**LECKEY**®

Code	Display Name	Safe Working Load
PED479240	Leckey Totstander, Ages 1-5 years	22 kg