

## Likoscale 350



- LikoScale is a convenient scale which is intended for weighing patients during lifting
- Together with Liko's lifts and accessories, LikoScale offers easy, comfortable and safe weighing
- After adjustment, the weight is shown in either kg or lbs on a display

 Contract approved

### Support @ Home

**Function:** Lifting persons

**Category:** Prescribed

Click [here](#) for more information

## Warranty Information

- 1 year

### Compliance

AS/NZS ISO 10535:2011



Code	Display Name	Safe Working Load
LSS393670	LikoScale 350	300kg
LSS393675	LikoScale 350 & Adaptor Kit	-
LSS393680	LikoScale Adapter Kit (contains Quick-release hook, Adapter 8mm and Q-link 13)	-