

TED - Anti-Embolic - Thigh Grip Top Stockings



- Thigh length with grip top
- Pressure break at knee, protecting the critical popliteal area and promoting free leg movement in the recumbent patient
- Available in 5 sizes, with choice of either, short, regular or long fit

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name	Colour
PTS645460	TED - Thigh Grip Top - Anti-embolic - Toe Hole - White	White
PTS641472	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - Lge / Short	White
PTS641473	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - XL / Short	White
PTS641474	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - XXL / Short	White
PTS641475	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - Med / Short	White
PTS641476	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - Sml / Reg	White
PTS641477	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - Med / Reg	White

Code	Display Name	Colour
PTS641478	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - Lge / Reg	White
PTS641479	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - XL / Reg	White
PTS641481	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - XXL / Reg	White
PTS641482	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - Sml / Long	White
PTS641483	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - Med / Long	White
PTS641485	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - XL / Long	White
PTS641486	TED Anti-embolic Stockings - Thigh w/ Grip top - Toe Hole - XXL / Long	White