

## **Large Capacity Medinizer**



**e-pill | Large Capacity Weekly Organiser and Medication Reminder** A practical and durable weekly medication management system with seven pillboxes and an integrated timer that allows up to six daily alarms. Designed for convenience at home, work, or while traveling — helping you stay on schedule and never miss a dose.



- Weekly Organiser: Includes seven individual pillboxes labeled Monday to Sunday for a full week of medication management.
- Four Daily Compartments: Each pillbox has four large sections labeled Morning, Noon, Evening, and Night to organise daily doses.
- Large Capacity: Each compartment holds up to 25 aspirin-sized tablets, ideal for multiple medications or supplements.
- Braille Markings: Each pillbox includes tactile labeling for accessibility and ease of use.
- **Durable Construction:** High-quality design with secure click-on lids and smooth sliding mechanism for easy access.
- **Built-in Timer:** Comes with a medication timer allowing up to six daily alarms to remind you when it's time to take your medication.
- **Easy Alarm Setup:** Simply slide the switches to set alarms on the hour; alarms automatically repeat daily.
- Audible and Visual Alerts: At alarm time, a 15-second beeping alarm sounds; if not stopped, it repeats after 5 minutes. A flashing indicator provides a visual alert for 5 minutes.
- **Simple to Use:** Press the stop button after taking your medication to silence the alarm.
- **Secure and Convenient:** Unique click-on lid design keeps pills safely stored while ensuring easy access when needed.
- Portable Design: Encased in a stylish vegan leather carry case for discretion



and protection.

Organize an entire week's medication with one setup: The e-pill Large
 Capacity Weekly Organizer simplifies complex medication schedules, ensures
 timely reminders, and offers portability with a discreet design.

## **Specifications:**

Canvas Bag Size: 6.5 in × 6.5 in × 2 in
Daily Pillbox Size: 3.5 in × 2 in × 0.75 in

• Total Weight: 15 oz

• Pill Capacity: 25 aspirin-sized tablets per compartment

• Color: Blue canvas bag with clear pillboxes

• Battery: 2 × AAA batteries (included)

• Battery Life: 6–12months, depending on the number of daily alarms set

Code Display Name

DLG275040 Large Capacity Medinizer Weekly Organizer Reminder