

Aspire Ultimate 8 Mattress Replacement



The Ultimate 8 mattress caters to users considered to be at high risk of pressure injury. Unique pump technology reduces noise, promoting improved sleep quality. The mattress is easy to clean and designed for durability.

- **HIGH LEVEL PRESSURE THERAPY:**
 - Promotes circulation and provides pressure redistribution to treat and prevent pressure injuries
 - Rapid release CPR allows quick deflation of the mattress, including head cells, in emergencies
 - Prevents bottoming out by increasing sacral zone air pressure when the mattress backrest is raised
 - 24 hr transport mode and static support layer create a security buffer allowing patients to remain on an inflated surface in the event of power failure
 - Easy-to-read LCD screen provides clear and precise information on cycle time, operating mode, and fault codes to assist with troubleshooting
 - Improves workplace safety via integrated power cord management in the mattress cover.
 - Designed for durability, the cell guard system improves mattress strength in key entry/exit areas
- **IMPROVES SLEEP QUALITY:**
 - Provides a whisper-quiet sleep environment using innovative pump to reduce noise and vibration
 - Improves sleep comfort by passing room temperature air beneath the sleep surface, helping to regulate body temperature
- **EASY INFECTION CONTROL:**

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

- Prevents fluid and bacterial ingress with welded seams and zip protection
- Multi-stretch cover reduces shear and offers a high moisture vapour transfer rate (MVTR) to control skin microclimate
- Easy to clean mattress cover can be laundered at up to 85°
- FEATURES OVERVIEW
 - Infection Control
 - Skin Microclimate
 - Alternating Pressure Relief
 - Easy Maintenance & Cleaning
 - Pressure Redistribution
 - Multiple Weight Settings

Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
PTM639320	Aspire Ultimate 8 Mattress Replacement - Single	2000 mm	900 mm	200 mm	30-180 kg	High Risk
PTM639330	Aspire Ultimate 8 Mattress Replacement - King Single	2000 mm	1050 mm	200 mm	30-180 kg	High Risk