

Jobst Bella Lite - Compression Hand - Ready To Wear



- Soft and air permeable
- Lightweight, knit construction
- Fits a large range of hand shapes
- Reduced compression at the wrist provides an effective and comfortable fit when the armsleeve and gauntlet are worn together
- Easy to put on, comfortable to wear
- Durable Applications:
 - Edema
 - Lymphoedema
 - Phlebitis
 - Post-thrombotic syndrome

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name
PTS642146	Compression Hand - Jobst Bella Lite - Ready To Wear - Gauntlet 20-30mmHg - Sz Sml
PTS642147	Compression Hand - Jobst Bella Lite - Ready To Wear - Gauntlet 20-30mmHg - Sz Med
PTS642148	Compression Hand - Jobst Bella Lite - Ready To Wear - Gauntlet 20-30mmHg - Sz Lge