

## Caeli Combination Mattress Replacement



- Caeli Combination 18 Air Cell Mattress Replacement with foam underlay
- Pillow head cells
- Fluid resistant cover with two-way stretch material
- Low air loss cells assist with cooling, reducing moisture build up.
- Quick pump connection
- Transport cap keeps mattress inflated for up to 24 hours
- 50mm foam underlay protects against bottoming out
- Cable management
- Dual CPR release valves for maximum efficiency in emergency situations
- Carry handles on both sides
- Safety strap allows mattress to be fixed to bed platform

✔ Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

## Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
PTM640100	Caeli Combination Mattress Replacement - 18 Cells - Single	2000 mm	840 mm	150 mm	180 kg	Medium Risk
PTM640110	Caeli Combination Mattress Replacement - 18 Cells - King Single	2000 mm	1080 mm	150 mm	180 kg	Medium Risk