

JAY - J3 Backrest



- The Jay-J3 series of backrests are highly customisable to suit your client's needs
- Customisations cover 5 spectrums: (1) Contour Depth; (2) Level of Support; (3) Back Support Height; (4) Backrest Width; (5) Mounting Hardware. Various options within each spectrum create an enormously flexible package
- We strongly recommend that you consult with your Aidacare C.A.S.E team member for any inquiries relating to these cushions
- The J3 Back is dedicated to eliminating compromises when different products are bundled. J3 ensures compatibility between solutions by encompassing a huge array of options within the one range
- The J3 Back features one of the widest ranges of available accessories. Multiple lateral sizes, Spine-Align positioning components, Whitmyer® headrests, Comfort Fit® chest straps, and trunk straps are an example of some of the available options
- Every J3 back features special new spacer fabric upholstery designed to wick moisture away from the skin and circulate air flow throughout the cover
- 135kg capacity on standard models and 225kg in models over 22"
- The J3 back provides up to 22° of back angle adjustment with minimal to zero loss of seat depth. This is especially helpful for individuals with postural deformities who may require more extreme angles
- Fitting the back to the user not the user to the back is the goal of the Spine-Align Posture Conforming system. This is accomplished by combining a wide range of individual foam components to create unique shapes that optimize postural support, evenly distribute pressure, and maximize function. Velcro strips attach to each shape so they can be placed into hundreds of configurations. The flexibility of the system is boundless

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Accessories for wheelchair seating

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_181003192_0105_1_2

Wheelchair Manual – Composite:

05_417_0105_1_2

Wheelchair Powered –

Composite: 05_419_0105_1_2

Click [here](#) for more information

Warranty Information

- 2 Years Backrest
- 6 Months Cover

Code	Display Name
JAY_J3_BACKREST	JAY - J3 Backrest
PTB533751	W`Chair Back - JAY J3 - Posterior & Lateral - Upper Thoracic - Med 18in
PTB533780	W`Chair Back - JAY J3 - Posterior & Lateral - Shoulder Height Small 18in
PTB533784	W`Chair Back - JAY J3 - Posterior & Deep Lateral - Shoulder Height Small 18in
PTB533838	W`Chair Back - JAY J3 - Laterals - Swing Away 6x6in