


Oxford Quickfit Deluxe Sling



- This sling is a development of the Quickfit design and is suitable for 95% of patients
- It can be used for amputees following assessment, and is available in sizes S to L
- FIM Scores: 1 and 2 and any other person from the floor

Warranty Information

- 1 Year



- ✓ Contracted DVA product
- ✓ Contract approved

ndis

Item Code:

05_123603164_0105_1_2

Sling – Composite:

05_421_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS371632	Quickfit Deluxe Sling - Polyester - Medium	227 kg	Polyester
LSS371700	Quikfit Deluxe Sling - Mesh - Medium	227 kg	Mesh