


## Oxford Quickfit Deluxe Sling



- This sling is a development of the Quickfit design and is suitable for 95% of patients
- It can be used for amputees following assessment, and is available in sizes S to L
- FIM Scores: 1 and 2 and any other person from the floor

### Warranty Information

- 1 Year

- ✔ Contracted DVA product
- ✔ Contract approved

### ndis

**Item Code:**

05\_123603164\_0105\_1\_2

**Sling – Composite:**

05\_421\_0105\_1\_2

Click [here](#) for more information



### Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS371632	Quickfit Deluxe Sling - Polyester - Medium	227 kg	Polyester
LSS371700	Quikfit Deluxe Sling - Mesh - Medium	227 kg	Mesh