


Sigvaris Magic Compression Stockings



- Sigvaris Magic are a versatile compression stocking
- Available in knee, thigh and waist high
- A delicate, slightly transparent stocking with a fine, elegant look
- Provides optimum leg support
- For every occasion and all ages
- Perfect fit and exceptional degree of comfort, thanks to high quality workmanship right down to the smallest of details
- Please specify fitting (open or closed toe) size and colour at time of ordering

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name	Colour
PTS649910	SIGVARIS Magic - Knee - 23-32mmHg XXS PLUS Norm Skin Open Toe	Skin
PTS649905	SIGVARIS Magic - Knee - 23-32mmHg XXS PLUS Norm Skin Closed Toe	Skin
PTS649900	SIGVARIS Magic - Knee - 23-32mmHg XXS PLUS Norm Black Open Toe	Black
PTS649885	SIGVARIS Magic - Knee - 23-32mmHg XXS PLUS Long Skin Closed Toe	Skin
PTS649890	SIGVARIS Magic - Knee - 23-32mmHg XXS PLUS Long Skin Open Toe	Skin
PTS649945	SIGVARIS Magic - Knee - 23-32mmHg L Norm Skin Closed Toe	Skin
PTS649880	SIGVARIS Magic - Knee - 23-32mmHg XXS PLUS Long Black Open Toe	Black

Code	Display Name	Colour
PTS649895	SIGVARIS Magic - Knee - 23-32mmHg XXS PLUS Norm Black Closed Toe	Black
PTS649915	SIGVARIS Magic - Knee - 23-32mmHg L Long Black Closed Toe	Black
PTS649920	SIGVARIS Magic - Knee - 23-32mmHg L Long Black Open Toe	Black
PTS649925	SIGVARIS Magic - Knee - 23-32mmHg L Long Skin Closed Toe	Skin
PTS649930	SIGVARIS Magic - Knee - 23-32mmHg L Long Skin Open Toe	Skin
PTS649875	SIGVARIS Magic - Knee - 23-32mmHg XXS PLUS Long Black Closed Toe	Black
PTS649940	SIGVARIS Magic - Knee - 23-32mmHg L Norm Black Open Toe	Black
PTS649830	SIGVARIS Magic - Knee - 23-32mmHg XL PLUS Norm Skin Open Toe	Skin
PTS649950	SIGVARIS Magic - Knee - 23-32mmHg L Norm Skin Open Toe	Skin
PTS649955	SIGVARIS Magic - Knee - 23-32mmHg M Long Black Closed Toe	Black
PTS649935	SIGVARIS Magic - Knee - 23-32mmHg L Norm Black Closed Toe	Black
PTS649820	SIGVARIS Magic - Knee - 23-32mmHg XL PLUS Norm Black Open Toe	Black
PTS649770	SIGVARIS Magic - Knee - 23-32mmHg S PLUS Long Skin Open Toe	Skin
PTS649775	SIGVARIS Magic - Knee - 23-32mmHg S PLUS Norm Black Closed Toe	Black
PTS649780	SIGVARIS Magic - Knee - 23-32mmHg S PLUS Norm Black Open Toe	Black
PTS649785	SIGVARIS Magic - Knee - 23-32mmHg S PLUS Norm Skin Closed Toe	Skin
PTS649790	SIGVARIS Magic - Knee - 23-32mmHg S PLUS Norm Skin Open Toe	Skin

Code	Display Name	Colour
PTS649795	SIGVARIS Magic - Knee - 23-32mmHg XL PLUS Long Black Closed Toe	Black
PTS649800	SIGVARIS Magic - Knee - 23-32mmHg XL PLUS Long Black Open Toe	Black
PTS649805	SIGVARIS Magic - Knee - 23-32mmHg XL PLUS Long Skin Closed Toe	Skin
PTS649840	SIGVARIS Magic - Knee - 23-32mmHg XS PLUS Long Black Open Toe	Black
PTS649815	SIGVARIS Magic - Knee - 23-32mmHg XL PLUS Norm Black Closed Toe	Black
PTS649870	SIGVARIS Magic - Knee - 23-32mmHg XS PLUS Norm Skin Open Toe	Skin
PTS649825	SIGVARIS Magic - Knee - 23-32mmHg XL PLUS Norm Skin Closed Toe	Skin
PTS649835	SIGVARIS Magic - Knee - 23-32mmHg XS PLUS Long Black Closed Toe	Black
PTS649960	SIGVARIS Magic - Knee - 23-32mmHg M Long Black Open Toe	Black
PTS649845	SIGVARIS Magic - Knee - 23-32mmHg XS PLUS Long Skin Closed Toe	Skin
PTS649850	SIGVARIS Magic - Knee - 23-32mmHg XS PLUS Long Skin Open Toe	Skin
PTS649855	SIGVARIS Magic - Knee - 23-32mmHg XS PLUS Norm Black Closed Toe	Black
PTS649860	SIGVARIS Magic - Knee - 23-32mmHg XS PLUS Norm Black Open Toe	Black
PTS649865	SIGVARIS Magic - Knee - 23-32mmHg XS PLUS Norm Skin Closed Toe	Skin
PTS649810	SIGVARIS Magic - Knee - 23-32mmHg XL PLUS Long Skin Open Toe	Skin
PTS650120	SIGVARIS Magic - Knee - 23-32mmHg XXS Long Black Open Toe	Black
PTS650000	SIGVARIS Magic - Knee - 23-32mmHg S Long Black Open Toe	Black
PTS650080	SIGVARIS Magic - Knee - 23-32mmHg XS Long Black Open Toe	Black

Code	Display Name	Colour
PTS650085	SIGVARIS Magic - Knee - 23-32mmHg XS Long Skin Closed Toe	Skin
PTS650090	SIGVARIS Magic - Knee - 23-32mmHg XS Long Skin Open Toe	Skin
PTS650095	SIGVARIS Magic - Knee - 23-32mmHg XS Norm Black Closed Toe	Black
PTS650100	SIGVARIS Magic - Knee - 23-32mmHg XS Norm Black Open Toe	Black
PTS650105	SIGVARIS Magic - Knee - 23-32mmHg XS Norm Skin Closed Toe	Skin
PTS650070	SIGVARIS Magic - Knee - 23-32mmHg XL Norm Skin Open Toe	Skin
PTS650115	SIGVARIS Magic - Knee - 23-32mmHg XXS Long Black Closed Toe	Black
PTS650065	SIGVARIS Magic - Knee - 23-32mmHg XL Norm Skin Closed Toe	Skin
PTS650125	SIGVARIS Magic - Knee - 23-32mmHg XXS Long Skin Closed Toe	Skin
PTS650130	SIGVARIS Magic - Knee - 23-32mmHg XXS Long Skin Open Toe	Skin
PTS650135	SIGVARIS Magic - Knee - 23-32mmHg XXS Norm Black Closed Toe	Black
PTS650140	SIGVARIS Magic - Knee - 23-32mmHg XXS Norm Black Open Toe	Black
PTS650145	SIGVARIS Magic - Knee - 23-32mmHg XXS Norm Skin Closed Toe	Skin
PTS650150	SIGVARIS Magic - Knee - 23-32mmHg XXS Norm Skin Open Toe	Skin
PTS649500	SIGVARIS Magic - Knee - 18-21mmHg M Norm Black Open Toe	Black
PTS649765	SIGVARIS Magic - Knee - 23-32mmHg S PLUS Long Skin Closed Toe	Skin
PTS650110	SIGVARIS Magic - Knee - 23-32mmHg XS Norm Skin Open Toe	Skin
PTS650020	SIGVARIS Magic - Knee - 23-32mmHg S Norm Black Open Toe	Black

Code	Display Name	Colour
PTS649970	SIGVARIS Magic - Knee - 23-32mmHg M Long Skin Open Toe	Skin
PTS649975	SIGVARIS Magic - Knee - 23-32mmHg M Norm Black Closed Toe	Black
PTS649980	SIGVARIS Magic - Knee - 23-32mmHg M Norm Black Open Toe	Black
PTS649985	SIGVARIS Magic - Knee - 23-32mmHg M Norm Skin Closed Toe	Skin
PTS649990	SIGVARIS Magic - Knee - 23-32mmHg M Norm Skin Open Toe	Skin
PTS649995	SIGVARIS Magic - Knee - 23-32mmHg S Long Black Closed Toe	Black
PTS650005	SIGVARIS Magic - Knee - 23-32mmHg S Long Skin Closed Toe	Skin
PTS650075	SIGVARIS Magic - Knee - 23-32mmHg XS Long Black Closed Toe	Black
PTS650015	SIGVARIS Magic - Knee - 23-32mmHg S Norm Black Closed Toe	Black
PTS649965	SIGVARIS Magic - Knee - 23-32mmHg M Long Skin Closed Toe	Skin
PTS650025	SIGVARIS Magic - Knee - 23-32mmHg S Norm Skin Closed Toe	Skin
PTS650030	SIGVARIS Magic - Knee - 23-32mmHg S Norm Skin Open Toe	Skin
PTS650035	SIGVARIS Magic - Knee - 23-32mmHg XL Long Black Closed Toe	Black
PTS650040	SIGVARIS Magic - Knee - 23-32mmHg XL Long Black Open Toe	Black
PTS650045	SIGVARIS Magic - Knee - 23-32mmHg XL Long Skin Closed Toe	Skin
PTS650050	SIGVARIS Magic - Knee - 23-32mmHg XL Long Skin Open Toe	Skin
PTS650055	SIGVARIS Magic - Knee - 23-32mmHg XL Norm Black Closed Toe	Black
PTS650060	SIGVARIS Magic - Knee - 23-32mmHg XL Norm Black Open Toe	Black

Code	Display Name	Colour
PTS650010	SIGVARIS Magic - Knee - 23-32mmHg S Long Skin Open Toe	Skin
PTS649400	SIGVARIS Magic - Knee - 18-21mmHg XXS PLUS Long Black Open Toe	Black
PTS649335	SIGVARIS Magic - Knee - 18-21mmHg XL PLUS Norm Black Closed Toe	Black
PTS649340	SIGVARIS Magic - Knee - 18-21mmHg XL PLUS Norm Black Open Toe	Black
PTS649345	SIGVARIS Magic - Knee - 18-21mmHg XL PLUS Norm Skin Closed Toe	Skin
PTS649350	SIGVARIS Magic - Knee - 18-21mmHg XL PLUS Norm Skin Open Toe	Skin
PTS649355	SIGVARIS Magic - Knee - 18-21mmHg XS PLUS Long Black Closed Toe	Black
PTS649360	SIGVARIS Magic - Knee - 18-21mmHg XS PLUS Long Black Open Toe	Black
PTS649365	SIGVARIS Magic - Knee - 18-21mmHg XS PLUS Long Skin Closed Toe	Skin
PTS649370	SIGVARIS Magic - Knee - 18-21mmHg XS PLUS Long Skin Open Toe	Skin
PTS649375	SIGVARIS Magic - Knee - 18-21mmHg XS PLUS Norm Black Closed Toe	Black
PTS649380	SIGVARIS Magic - Knee - 18-21mmHg XS PLUS Norm Black Open Toe	Black
PTS649385	SIGVARIS Magic - Knee - 18-21mmHg XS PLUS Norm Skin Closed Toe	Skin
PTS649760	SIGVARIS Magic - Knee - 23-32mmHg S PLUS Long Black Open Toe	Black
PTS649395	SIGVARIS Magic - Knee - 18-21mmHg XXS PLUS Long Black Closed Toe	Black
PTS649320	SIGVARIS Magic - Knee - 18-21mmHg XL PLUS Long Black Open Toe	Black
PTS649405	SIGVARIS Magic - Knee - 18-21mmHg XXS PLUS Long Skin Closed Toe	Skin
PTS649410	SIGVARIS Magic - Knee - 18-21mmHg XXS PLUS Long Skin Open Toe	Skin

Code	Display Name	Colour
PTS649415	SIGVARIS Magic - Knee - 18-21mmHg XXS PLUS Norm Black Closed Toe	Black
PTS649420	SIGVARIS Magic - Knee - 18-21mmHg XXS PLUS Norm Black Open Toe	Black
PTS649425	SIGVARIS Magic - Knee - 18-21mmHg XXS PLUS Norm Skin Closed Toe	Skin
PTS649430	SIGVARIS Magic - Knee - 18-21mmHg XXS PLUS Norm Skin Open Toe	Skin
PTS649435	SIGVARIS Magic - Knee - 18-21mmHg L Long Black Closed Toe	Black
PTS649440	SIGVARIS Magic - Knee - 18-21mmHg L Long Black Open Toe	Black
PTS649445	SIGVARIS Magic - Knee - 18-21mmHg L Long Skin Closed Toe	Skin
PTS649450	SIGVARIS Magic - Knee - 18-21mmHg L Long Skin Open Toe	Skin
PTS649510	SIGVARIS Magic - Knee - 18-21mmHg M Norm Skin Open Toe	Skin
PTS649460	SIGVARIS Magic - Knee - 18-21mmHg L Norm Black Open Toe	Black
PTS649390	SIGVARIS Magic - Knee - 18-21mmHg XS PLUS Norm Skin Open Toe	Skin
PTS649260	SIGVARIS Magic - Knee - 18-21mmHg M PLUS Norm Black Open Toe	Black
PTS649195	SIGVARIS Magic - Knee - 18-21mmHg L PLUS Long Black Closed Toe	Black
PTS649200	SIGVARIS Magic - Knee - 18-21mmHg L PLUS Long Black Open Toe	Black
PTS649205	SIGVARIS Magic - Knee - 18-21mmHg L PLUS Long Skin Closed Toe	Skin
PTS649210	SIGVARIS Magic - Knee - 18-21mmHg L PLUS Long Skin Open Toe	Skin
PTS649215	SIGVARIS Magic - Knee - 18-21mmHg L PLUS Norm Black Closed Toe	Black
PTS649220	SIGVARIS Magic - Knee - 18-21mmHg L PLUS Norm Black Open Toe	Black

Code	Display Name	Colour
PTS649225	SIGVARIS Magic - Knee - 18-21mmHg L PLUS Norm Skin Closed Toe	Skin
PTS649230	SIGVARIS Magic - Knee - 18-21mmHg L PLUS Norm Skin Open Toe	Skin
PTS649235	SIGVARIS Magic - Knee - 18-21mmHg M PLUS Long Black Closed Toe	Black
PTS649240	SIGVARIS Magic - Knee - 18-21mmHg M PLUS Long Black Open Toe	Black
PTS649245	SIGVARIS Magic - Knee - 18-21mmHg M PLUS Long Skin Closed Toe	Skin
PTS649330	SIGVARIS Magic - Knee - 18-21mmHg XL PLUS Long Skin Open Toe	Skin
PTS649255	SIGVARIS Magic - Knee - 18-21mmHg M PLUS Norm Black Closed Toe	Black
PTS649325	SIGVARIS Magic - Knee - 18-21mmHg XL PLUS Long Skin Closed Toe	Skin
PTS649265	SIGVARIS Magic - Knee - 18-21mmHg M PLUS Norm Skin Closed Toe	Skin
PTS649270	SIGVARIS Magic - Knee - 18-21mmHg M PLUS Norm Skin Open Toe	Skin
PTS649275	SIGVARIS Magic - Knee - 18-21mmHg S PLUS Long Black Closed Toe	Black
PTS649280	SIGVARIS Magic - Knee - 18-21mmHg S PLUS Long Black Open Toe	Black
PTS649285	SIGVARIS Magic - Knee - 18-21mmHg S PLUS Long Skin Closed Toe	Skin
PTS649290	SIGVARIS Magic - Knee - 18-21mmHg S PLUS Long Skin Open Toe	Skin
PTS649295	SIGVARIS Magic - Knee - 18-21mmHg S PLUS Norm Black Closed Toe	Black
PTS649300	SIGVARIS Magic - Knee - 18-21mmHg S PLUS Norm Black Open Toe	Black
PTS649305	SIGVARIS Magic - Knee - 18-21mmHg S PLUS Norm Skin Closed Toe	Skin
PTS649310	SIGVARIS Magic - Knee - 18-21mmHg S PLUS Norm Skin Open Toe	Skin

Code	Display Name	Colour
PTS649315	SIGVARIS Magic - Knee - 18-21mmHg XL PLUS Long Black Closed Toe	Black
PTS649465	SIGVARIS Magic - Knee - 18-21mmHg L Norm Skin Closed Toe	Skin
PTS649250	SIGVARIS Magic - Knee - 18-21mmHg M PLUS Long Skin Open Toe	Skin
PTS649695	SIGVARIS Magic - Knee - 23-32mmHg L PLUS Norm Black Closed Toe	Black
PTS649630	SIGVARIS Magic - Knee - 18-21mmHg XS Norm Skin Open Toe	Skin
PTS649455	SIGVARIS Magic - Knee - 18-21mmHg L Norm Black Closed Toe	Black
PTS649640	SIGVARIS Magic - Knee - 18-21mmHg XXS Long Black Open Toe	Black
PTS649645	SIGVARIS Magic - Knee - 18-21mmHg XXS Long Skin Closed Toe	Skin
PTS649650	SIGVARIS Magic - Knee - 18-21mmHg XXS Long Skin Open Toe	Skin
PTS649655	SIGVARIS Magic - Knee - 18-21mmHg XXS Norm Black Closed Toe	Black
PTS649660	SIGVARIS Magic - Knee - 18-21mmHg XXS Norm Black Open Toe	Black
PTS649665	SIGVARIS Magic - Knee - 18-21mmHg XXS Norm Skin Closed Toe	Skin
PTS649670	SIGVARIS Magic - Knee - 18-21mmHg XXS Norm Skin Open Toe	Skin
PTS649675	SIGVARIS Magic - Knee - 23-32mmHg L PLUS Long Black Closed Toe	Black
PTS649680	SIGVARIS Magic - Knee - 23-32mmHg L PLUS Long Black Open Toe	Black
PTS649470	SIGVARIS Magic - Knee - 18-21mmHg L Norm Skin Open Toe	Skin
PTS649690	SIGVARIS Magic - Knee - 23-32mmHg L PLUS Long Skin Open Toe	Skin
PTS649625	SIGVARIS Magic - Knee - 18-21mmHg XS Norm Skin Closed Toe	Skin

Code	Display Name	Colour
PTS649700	SIGVARIS Magic - Knee - 23-32mmHg L PLUS Norm Black Open Toe	Black
PTS649705	SIGVARIS Magic - Knee - 23-32mmHg L PLUS Norm Skin Closed Toe	Skin
PTS649710	SIGVARIS Magic - Knee - 23-32mmHg L PLUS Norm Skin Open Toe	Skin
PTS649715	SIGVARIS Magic - Knee - 23-32mmHg M PLUS Long Black Closed Toe	Black
PTS649720	SIGVARIS Magic - Knee - 23-32mmHg M PLUS Long Black Open Toe	Black
PTS649725	SIGVARIS Magic - Knee - 23-32mmHg M PLUS Long Skin Closed Toe	Skin
PTS649730	SIGVARIS Magic - Knee - 23-32mmHg M PLUS Long Skin Open Toe	Skin
PTS649735	SIGVARIS Magic - Knee - 23-32mmHg M PLUS Norm Black Closed Toe	Black
PTS649740	SIGVARIS Magic - Knee - 23-32mmHg M PLUS Norm Black Open Toe	Black
PTS649745	SIGVARIS Magic - Knee - 23-32mmHg M PLUS Norm Skin Closed Toe	Skin
PTS649750	SIGVARIS Magic - Knee - 23-32mmHg M PLUS Norm Skin Open Toe	Skin
PTS649755	SIGVARIS Magic - Knee - 23-32mmHg S PLUS Long Black Closed Toe	Black
PTS649685	SIGVARIS Magic - Knee - 23-32mmHg L PLUS Long Skin Closed Toe	Skin
PTS649545	SIGVARIS Magic - Knee - 18-21mmHg S Norm Skin Closed Toe	Skin
PTS649475	SIGVARIS Magic - Knee - 18-21mmHg M Long Black Closed Toe	Black
PTS649480	SIGVARIS Magic - Knee - 18-21mmHg M Long Black Open Toe	Black
PTS649485	SIGVARIS Magic - Knee - 18-21mmHg M Long Skin Closed Toe	Skin
PTS649490	SIGVARIS Magic - Knee - 18-21mmHg M Long Skin Open Toe	Skin

Code	Display Name	Colour
PTS649495	SIGVARIS Magic - Knee - 18-21mmHg M Norm Black Closed Toe	Black
PTS649505	SIGVARIS Magic - Knee - 18-21mmHg M Norm Skin Closed Toe	Skin
PTS649515	SIGVARIS Magic - Knee - 18-21mmHg S Long Black Closed Toe	Black
PTS649520	SIGVARIS Magic - Knee - 18-21mmHg S Long Black Open Toe	Black
PTS649525	SIGVARIS Magic - Knee - 18-21mmHg S Long Skin Closed Toe	Skin
PTS649530	SIGVARIS Magic - Knee - 18-21mmHg S Long Skin Open Toe	Skin
PTS649635	SIGVARIS Magic - Knee - 18-21mmHg XXS Long Black Closed Toe	Black
PTS649540	SIGVARIS Magic - Knee - 18-21mmHg S Norm Black Open Toe	Black
PTS649620	SIGVARIS Magic - Knee - 18-21mmHg XS Norm Black Open Toe	Black
PTS649590	SIGVARIS Magic - Knee - 18-21mmHg XL Norm Skin Open Toe	Skin
PTS649615	SIGVARIS Magic - Knee - 18-21mmHg XS Norm Black Closed Toe	Black
PTS649610	SIGVARIS Magic - Knee - 18-21mmHg XS Long Skin Open Toe	Skin
PTS649605	SIGVARIS Magic - Knee - 18-21mmHg XS Long Skin Closed Toe	Skin
PTS649535	SIGVARIS Magic - Knee - 18-21mmHg S Norm Black Closed Toe	Black
PTS649595	SIGVARIS Magic - Knee - 18-21mmHg XS Long Black Closed Toe	Black
PTS649550	SIGVARIS Magic - Knee - 18-21mmHg S Norm Skin Open Toe	Skin
PTS649585	SIGVARIS Magic - Knee - 18-21mmHg XL Norm Skin Closed Toe	Skin
PTS649580	SIGVARIS Magic - Knee - 18-21mmHg XL Norm Black Open Toe	Black

Code	Display Name	Colour
PTS649575	SIGVARIS Magic - Knee - 18-21mmHg XL Norm Black Closed Toe	Black
PTS649570	SIGVARIS Magic - Knee - 18-21mmHg XL Long Skin Open Toe	Skin
PTS649565	SIGVARIS Magic - Knee - 18-21mmHg XL Long Skin Closed Toe	Skin
PTS649560	SIGVARIS Magic - Knee - 18-21mmHg XL Long Black Open Toe	Black
PTS649555	SIGVARIS Magic - Knee - 18-21mmHg XL Long Black Closed Toe	Black
PTS649600	SIGVARIS Magic - Knee - 18-21mmHg XS Long Black Open Toe	Black