

Liko UltraSling (Mod. 06)



- Designed for safe and comfortable lifting of bariatric patients
- Supports the entire back up to the shoulders
- Leg supports are installed individually on the sling bar, which gives the stomach space and the legs optimal support
- Available in Short, Regular and Long lengths
- Available in Mesh & Polyester

Warranty Information

- 1 Year



 Contract approved

ndis

Item Code:

05_123603166_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394410	Ultrasling Mod 6 Sling - Reinforced Leg Support - Polyester - REGULAR	500 kg	Polyester (Loops)