

Jobst Body Adhesive



- Helps Jobst garments stay in place. Ideal for thigh length stockings, when garter belt or silicone bands are undesirable
- Easy to apply and wash off
- Can be used with crew socks, knee high socks and stockings, thigh high stockings, any compression garment, even shoulder straps, wigs, toupees, and other orthopaedic & surgical devices

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information



Code	Display Name
PTS641100	Jobst - Body Adhesive - 59ml Plastic Bottle