

Aspire Raphael X CUSTOM Quattro Lift Recline Chair

Customise the Aspire Raphael X Lift Chair to match your unique style and house
Available in a range of commercial quality fabric for long lasting durability



HEALTHCARE VINYL LUSTRELL RANGER RANGE



- **Soft Touch Fabric** - The soft touch weave and diverse colour palette makes the Manisa range a highly practical and functional solution. These fabrics are comforting, easy to clean and durable, yet still soft to touch.
- **Leather** - A luxurious artisan leather with a soft glove-like texture, will retain its unique suppleness, provided they are regularly cared for. The leather is produced and tested to international standards for properties such as rubbing, flex resistance, light fastness, thickness, and strength.
- **Microsuede** - A Suede like texture combined with high performance characteristics ensures the Eastwood fabric range remains easy to clean and maintain.
- **Healthcare Vinyl** The Lustrell Ranger range is high performance vinyl with HealthGuard Antibacterial and Antimicrobial treatment to assist with infection control and hygiene.
- Custom order lead time varies, approximately 5 - 6 weeks depending on fabric availability.
- Custom orders are subject to terms and conditions, please see page two on the order form below.
- All colour swatches shown may not be a true representation of the actual product and should be used as a guide only.
- All colour swatches shown may not be a true representation of the actual product and should be used as a guide only.

Raphael X lift recline chair details

Designed for enhanced comfort and customisable support. The Raphael X delivers exceptional comfort with features to suit a diverse range of users. The adjustable backrest includes access for personalised postural adjustment and lateral support,

providing improved upright sitting. Its plush design compliments any modern home environment.

- **Built-in Lateral Supports:** Helps maintain an upright sitting posture, improving comfort and stability
- **Integrated Footrest Extension:** Helps with fatigue and improves lower limb comfort*
- **Adjustable Seat Depth:** Offers 25mm of depth adjustment, allowing taller users to achieve a better-seated position
- **Protective Covers for High-Wear Areas:** Removable covers for the headrest and armrests to keep your chair looking new
- **Refined, Contemporary Design:** Larger, plush armrests provide extra comfort while maintaining a streamlined look that suits home environment
- **All-day Battery Backup:** A lithium-ion battery ensures continued operation during power outages, supporting safety and independence
- **Classic Backrest Option:** Available as an additional accessory for users who prefer a more traditional support profile *The Petite model is designed without a leg extension, catering to the needs of shorter users.

Warranty Information

- 2 Years

Code	Display Name	RRP	Box Quantity	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
CHP227640FA	Raphael X CUSTOM Quattro Lift Recline Chair - PETITE - Soft Touch Fabric	4159	1	820 mm	1060 mm	130 kg	450 - 475 mm	470 mm	450 mm	450 mm

Code	Display Name	RRP	Box Quantity	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
CHP227640LT	Raphael X CUSTOM Quattro Lift Recline Chair - PETITE - Leather	4779	1	820 mm	1060 mm	130 kg	450 - 475 mm	470 mm	450 mm	450 mm
CHP227640MS	Raphael X CUSTOM Quattro Lift Recline Chair - PETITE - Microsuede	4159	1	820 mm	1060 mm	130 kg	450 - 475 mm	470 mm	450 mm	450 mm
CHP227640VL	Raphael X CUSTOM Quattro Lift Recline Chair - PETITE - Healthcare Vinyl	4159	1	820 mm	1060 mm	130 kg	450 - 475 mm	470 mm	450 mm	450 mm
CHP227650FA	Raphael X CUSTOM Quattro Lift Recline Chair - SMALL Soft Touch Fabric	4159	1	850 mm	1060 mm	130 kg	480 - 505 mm	500 mm	450 mm	450 mm
CHP227650LT	Raphael X CUSTOM Quattro Lift Recline Chair - SMALL - Leather	4779	1	850 mm	1060 mm	130 kg	480 - 505 mm	500 mm	450 mm	450 mm
CHP227650MS	Raphael X CUSTOM Quattro Lift Recline Chair - SMALL - Microsuede	4159	1	850 mm	1060 mm	130 kg	480 - 505 mm	500 mm	450 mm	450 mm
CHP227650VL	Raphael X CUSTOM Quattro Lift Recline Chair - SMALL - Healthcare Vinyl	4159	1	850 mm	1060 mm	130 kg	480 - 505 mm	500 mm	450 mm	450 mm
CHP227660FA	Raphael X CUSTOM Quattro Lift Recline Chair - MEDIUM Soft Touch Fabric	4159	1	875 mm	1100 mm	130 kg	500 - 525 mm	525 mm	470 mm	470 mm
CHP227660LT	Raphael X CUSTOM Quattro Lift Recline Chair - MEDIUM - Leather	4779	1	875 mm	1100 mm	130 kg	500 - 525 mm	525 mm	470 mm	470 mm
CHP227660MS	Raphael X CUSTOM Quattro Lift Recline Chair - MEDIUM - Microsuede	4159	1	875 mm	1100 mm	130 kg	500 - 525 mm	525 mm	470 mm	470 mm
CHP227660VL	Raphael X CUSTOM Quattro Lift Recline Chair - MEDIUM - Healthcare Vinyl	4159	1	875 mm	1100 mm	130 kg	500 - 525 mm	525 mm	470 mm	470 mm

Code	Display Name	RRP	Box Quantity	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
CHP227670FA	Raphael X CUSTOM Quattro Lift Recline Chair - LARGE Soft Touch Fabric	4159	1	900 mm	1100 mm	130 kg	530 - 555 mm	550 mm	500 mm	500 mm
CHP227670LT	Raphael X CUSTOM Quattro Lift Recline Chair - LARGE - Leather	4779	1	900 mm	1100 mm	130 kg	530 - 555 mm	550 mm	500 mm	500 mm
CHP227670MS	Raphael X CUSTOM Quattro Lift Recline Chair - LARGE - Microsuede	4159	1	900 mm	1100 mm	130 kg	530 - 555 mm	550 mm	500 mm	500 mm
CHP227670VL	Raphael X CUSTOM Quattro Lift Recline Chair - LARGE - Healthcare Vinyl	4159	1	900 mm	1100 mm	130 kg	530 - 555 mm	550 mm	500 mm	500 mm
CHP227680FA	Raphael X CUSTOM Quattro Lift Recline Chair - MAXI- Soft Touch Fabric	4369	1	1020 mm	1120 mm	180 kg	540 - 565 mm	650 mm	510 mm	510 mm
CHP227680LT	Raphael X CUSTOM Quattro Lift Recline Chair - MAXI - Leather	4989	1	1020 mm	1120 mm	180 kg	540 - 565 mm	650 mm	510 mm	510 mm
CHP227680MS	Raphael X CUSTOM Quattro Lift Recline Chair - MAXI - Microsuede	4369	1	1020 mm	1120 mm	180 kg	540 - 565 mm	650 mm	510 mm	510 mm
CHP227680VL	Raphael X CUSTOM Quattro Lift Recline Chair - MAXI - Healthcare Vinyl	4369	1	1020 mm	1120 mm	180 kg	540 - 565 mm	650 mm	510 mm	510 mm