

OPPO 2031 HINGED KNEE OPEN PATELLA SUPPORT



- May be used for stabilization and support of the knee following strains, injury or surgery
- Two hook and loop straps insure comfortable fit
- Elastic brace with motion hinges provides medial and lateral support for the unstable knee
- For arthritic knee conditions, minor sprains and strains
- Provides firm support for the maintenance of joint structures with gentle compression over the kneecap Application:
- Weak knee
- Arthritis
- Knee instability

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS308370	OPPO 2031 HINGED KNEE OPEN PATELLA SUPPORT - XXL