

eCart Accessories



Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width
WAT901115	eCart - Drive Unit	-	-
WAT902110	Basket	-	-
WAT902090	Accessory Rail - for eCart	-	-
WAT902120	Bottle Cabinet - Fitted with Key Lock	-	-
WAT902130	Cardiac Board	-	-
WAT902140	Chart Holder	-	-
WAT902150	Dividers (1 x Set) (4 x 3 = 12)	-	-

Code	Display Name	Total Length	Total Width
WAT902160	Divider Name Tabs	-	-
WAT902170	Defibrillator Arm	-	-
WAT902220	Equipment Shelf - Double	-	-
WAT902230	Equipment Shelf - Single	-	-
WAT902240	Folding Side Shelf	-	-
WAT902250	Guard Rail	-	-
WAT902260	Intravenous Pole with Corner Bracket	-	-
WAT902280	Oxygen Bottle Carrier	-	-
WAT902300	Sharps Holder	-	-
WAT902310	Stainless Steel Top	-	-
WAT902320	Tablet Arm	-	-
WAT902330	Tilt Container - 4	-	-
WAT902340	Tilt container - 5	-	-
WAT902350	Tilt container - 6	-	-
WAT902360	Tilt container - 9	-	-
WAT902370	Mount Arms for Tilt Container (Pair)	-	-
WAT902380	Waste Bin	-	-

Code	Display Name	Total Length	Total Width
WAT902390	Webster Clips (Pair) to Suit Wire Webster Rack	-	-
WAT902400	Glove Dispenser	-	-
WAT902105	Anti Tamper Lock	-	-
WAT902410	Buffer corner for eCart	100 mm	-
WAT902420	Drawer runner for eCart range	402 mm	45 mm