

Conni Reusable Bed Pad



Conni.

- The Conni absorbent and waterproof bed pads assist in the management of incontinence, bed wetting, and excessive perspiration
- Designed to hold 2 litres over an eight hour period, which means they are the perfect addition to any bed that requires full void protection
- Cost-effective and environmentally sustainable alternative to disposable
- Certified as being free from harmful materials and are safe even against newborn skin
- Soft and breathable
- Washable and reusable

Warranty Information

- 1 Year



ndis

Item Code:

03_092489060_0103_1_1

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Colour	Capacity
CON228520	Conni Bed Pad - Mauve	63.5	11	950 mm	850 mm	Mauve	2 litres
CON228525	Conni Bed Pad - Teal Blue	63.5	11	950 mm	850 mm	Teal Blue	2 litres
CON228526	Conni Bed Pad - White	63.5	-	950 mm	850 mm	White	2 litres
CON228528	Conni Bed Pad - Animal print	63.5	-	950 mm	850 mm	Animal print	2 litres
CON228518	Conni Bed Pad - Charcoal	63.5	-	850 mm	950 mm	Charcoal	2 litres