

## Aidacare Aspire A200S Aluminium Standing Lifter



The Aidacare Aspire A200S is a compact, lightweight and manoeuvrable standing lifter ideal for home, community and facility environments.

- **Adjustable Depth Yoke:** Yoke is depth adjustable for optimal patient fit
- **Ergonomic Handles:** Provide multiple grip points to suit carers of different heights and strengths
- **Ergonomic Handset:** The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- **Powered Leg Function:** Opens and closes the legs for lifting stability or accessibility
- **Build Quality:** Strong powder coated aluminium frame with industry leading Linak electronics and warranty provision
- **Safety:** Critical features include emergency stop button, emergency lowering and antientrapment technology, with both patient and caregiver safety at heart
- **Adjustable Padded Legrest:** Provide comfort and support during transfers & can be adjusted to suit the user
- **Intelligent Diagnostics:** State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- **Tente Casters:** 2 x 75 mm front castors for maximum under bed clearance and accessibility. 2 x 100 mm rear locking castors for safety and manoeuvrability.
- Item Weight: 30.4 kg
- Tested and compliant to AS/NZS ISO 10535:2011 (the Australian standard for hoists)

✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Lifting persons

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_123604167\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS ISO 10535:2011

## Warranty Information

- Frame: 5 Years
- Electronics: 3 Years
- Battery and Handset: 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Electric Leg	Safe Working Load	Material
LSS390640	Aidacare Aspire A200S Standing Patient Lifter	985 mm	590 - 915 mm	1455 mm	890 - 1395 mm	490 - 815 mm	True	200 kg	Aluminium