

Aidacare Aspire 200S Steel Standing Lifter



The Aidacare Aspire 200S is a robust and durable steel standing lifter capable in the most demanding facility environments

- **Bullhorn Lifting** Handles for patient stability and 6 sling attachment points
- **Ergonomic Handles** Provide multiple grip points to suit carers of different heights and strengths
- **Ergonomic Handset** The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- **Powered Leg Function** Opens and closes the legs for lifting stability or accessibility
- **Sturdy steel frame** with industry leading Linak electronics and warranty provision
- **Safety** Critical features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart
- **Adjustable Padded Legrests** Provide comfort and support during transfers & can be adjusted to suit the patient
- **Intelligent Diagnostics** State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- **Tente Casters:** 2 x 75 mm front castors for maximum under bed clearance and accessibility. 2 x 100 mm rear locking castors for safety and manoeuvrability
- Item Weight: 44.6 kg
- Tested and compliant to AS/NZS ISO 10535:2011 (the Australian standard for hoists)

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123604167_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Warranty Information

- Frame: 5 Years
- Electronics: 3 Years
- Battery and Handset: 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Lift Range	Leg Gait Range	Electric Leg
LSS390650	Aidacare Aspire 200S Standing Patient Lifter	1150 mm	656 - 1016 mm	1790 mm	200 kg	Steel	44.6 kg	850 - 1740 mm	575 - 936 mm	True