

The Procedure Cart



- Adjustable push handle, quick and easy to enable users to set push handle at the correct ergonomic height
- Storage trays
- Adjustable intravenous pole
- Recessed ABS formed plastic top
- Slide out shelf fitted with runners for smooth stable operation
- Storage basket
- Waste bin
- Sharps holder
- Corner buffers
- Diameter 125 anti-static total lock castors
- Large bottom drawer
- Raised drawer handles designed for grip and easy cleaning
- Epoxy powder coated
- 5 drawers fitted with soft close runners and damper prevents drawers from opening while in motion



Warranty Information

- 5 Years Frame

Code	Display Name	Total Length	Total Width	Total Height
WAT901220	Procedure Cart (Grey) - 5 Drawer + IV Pole + Storage Basket	520 mm	760 mm	1045 mm