

## Wonder Sheet - Fitted



- The 'Wonder Sheet' is an effective therapeutic tool because it can make your mobility in bed easier.
- The satin panel, which is 90cm long, sits under your shoulders, back and hips and allows you to slide your hips over more easily without having to lift them.
- Helpful when getting in and out of bed because you can pivot on your hips without worrying about the friction from normal cotton sheets.
- This unique design creatively combines 'Glossy Charmeuse' satin and a 50/50 blend cotton/polyester fabric which is beautiful to sleep on while being very durable.
- The 350 mm depth is suitable for mattresses up to 300 mm in height.

 Contract approved

### Support @ Home

**Function:** Beds and bed equipment

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

03\_092100082\_0103\_1\_1

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height
BEA019940	Wonder Sheet - Fitted - Single - 35cm Depth	1930 mm	910 mm	350 mm
BEA019950	Wonder Sheet - Fitted - Single - 45cm Depth	1930 mm	910 mm	450 mm
BEA019920	Wonder Sheet - Fitted - Long Single - 35cm Depth	2030 mm	910 mm	350 mm
BEA019890	Wonder Sheet - Fitted - King Single - 35cm Depth	2030 mm	1070 mm	350 mm
BEA019932	Wonder Sheet - Fitted - King Single - 45cm Depth	2030 mm	1070 mm	450 mm
BEA019860	Wonder Sheet - Fitted - Double - 35cm Depth	1930 mm	1350 mm	350 mm
BEA019930	Wonder Sheet - Fitted - Queen - 35cm Depth	2030 mm	1800 mm	350 mm

Code	Display Name	Total Length	Total Width	Total Height
BEA019931	Wonder Sheet - Fitted - Queen - 45cm Depth	2030 mm	1520 mm	450 mm
BEA019933	Wonder Sheet - Fitted - King - 45cm Depth	2030 mm	1800 mm	450 mm
BEA019960	Wonder Sheet - Pillow Case - 1/2 and 1/2 (pair)	700 mm	450 mm	-