

## NuboMed X5000 Medication Dispensing Cart



NuboMed

- **Material:** Constructed from anti-microbial ABS and aluminium alloy
- **Integrated Design:** The handle and keyboard tray are seamlessly integrated from a single piece of antimicrobial ABS material, ensuring a gap-free structure with a spill-proof groove design
- **Double-Pillar Design:** Drawer housing securely mounted to dual aluminium alloy pillars
- **Expanded Work Surface:** Pulls out from the rear handle for convenient use
- **Monitor Mount:** Height adjustable from 0-150mm, with 90° left and right swivel, 30° upward tilt, and portrait-to-landscape rotation
- **Button:** PC power on/off button, screen backlight on/off switch
- **Drawer Housing:** ABS drawer housing with lock, includes four configurable aluminium alloy drawers (1U+2U) with internal dividers
- **Medication Compartment:** Equipped with 2 sets of medication compartments, each containing 6 resin medication boxes made through a one-piece moulding process. They support electronic lock control and automatic ejection, and also feature an emergency unlocking function
- **Base:** Designed to prevent tripping and made from a single-piece of aluminium alloy with an antibacterial powder coating finish
- **Castors:** Equipped with 2x 4" lockable and 2x 4" tracking castors
- **Ideal For:** Closed Loop Intelligent Medication Management



### Warranty Information

- 2 Years

Code	Display Name
WAT908170	NuboMed - X5000 - Medication Dispensing Cart