

2 Section Fixed Examination Couch



- Fixed height 82 cm
- New easy lift and lock mechanism for easy assembly
- Lift up back section for assessment
- Comes flat packed to save on freight
- Padding thickness 5 cm
- Foam density: 36g/m2
- Breathe hole available on request



Compliance

IEC 60601-2-52:2009

Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Colour
WAT821575	Exam Couch - 2 Section Fixed Manual - Black	1950 mm	710 mm	820 mm	280 kg	Steel/Vinyl/Foam	Black