

Malvestio ONEday Electric Stretcher



- 4-section fully electric stretcher offering the functions of a bed with the smaller footprint and mobility of a stretcher
- Electric trendelenburg and reverse-trendelenburg functions to meet clinical and patient handling needs
- 200 mm castors minimise rolling resistance
- Fifth central wheel improves steering and reduces manual handling burden in transportation
- 260 kg SWL and inbuilt leg extension facilitate changing patient and clinical needs across multiple ward requirements
- Hand Control flexible arm promotes patient independence and comfort
- Enhanced battery-pack, with day-long operating time, dedicated battery charger for a full and fast recharge
- Audible alarm system sounds when mains power is lost and when brakes are unlocked
- Quick Release CPR for back rest and knee break ensures the stretcher is emergency ready at all times
- Frame is finished with Biocote, a technologically advanced surface coating that promotes infection control via its anti-bacterial properties
- Back rest angle 0 to 70 degrees and knee bend angle 0 to 30 degrees

✓ Contract approved

Compliance

IEC 60601-2-52:2009

Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
TSE674100	Malvestio ONEday Electric Hospital Stretcher	1990 - 2170 mm	890 mm	400 - 830 mm	260 kg