

## Venosan 5000 Compression Stockings



Venosan 5000 Compression Stockings are the premium choice for those seeking superior quality. Featuring the magic of silver with pure silver-coated X-Static fiber, these stockings maintain silver's properties even after 150 wash cycles. The silver's excellent heat conductivity, combined with microfiber TACTEL, rapidly dissipates skin moisture, making them durable and ideal for active individuals.

- Pure silver-coated X-Static fiber retains properties after 150 wash cycles.
- Silver's heat conductive properties in combination with TACTEL microfiber dissipate skin moisture.
- Durable construction with an innovative yarn blend of silver and cotton, suitable for active men and women.
- **COMPOSITION:** 40% TACTEL Polyamide, 30% Elastane, 20% Cotton, 10% X-Static - The Silver Fiber

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Supporting and stimulating blood circulation

**Category:** Low risk

Click [here](#) for more information

Code	Display Name	Colour
FBS308690	Venosan 5001 (Nylon) - Knee - 18-22MmHg - Unisex - O/Toe - SML - Long - Mexico	Mexico
FBS308700	Venosan 5001 (Nylon) - Knee - 18-22MmHg - Unisex - O/Toe - MED - Long - Mexico	Mexico
FBS308710	Venosan 5001 (Nylon) - Knee - 18-22MmHg - Unisex - O/Toe - LGE - Long - Mexico	Mexico
FBS308720	Venosan 5001 (Nylon) - Knee - 18-22MmHg - Unisex - O/Toe - XLGE - Long - Mexico	Mexico
FBS308730	Venosan 5001 (Nylon) - Knee - 18-22MmHg - Unisex - C/Toe - SML - Long - Mexico	Mexico
FBS308740	Venosan 5001 (Nylon) - Knee - 18-22MmHg - Unisex - C/Toe - MED - Long - Mexico	Mexico

Code	Display Name	Colour
FBS308750	Venosan 5001 (Nylon) - Knee - 18-22Mmhg - Unisex - C/Toe - LGE - Long - Mexico	Mexico
FBS308760	Venosan 5001 (Nylon) - Knee - 18-22Mmhg - Unisex - C/Toe - XLGE - Long - Mexico	Mexico
FBS308770	Venosan 5002 (Nylon) - Knee - 23-32Mmhg - Unisex - O/Toe - SML - Long - Mexico	Mexico
FBS308780	Venosan 5002 (Nylon) - Knee - 23-32Mmhg - Unisex - O/Toe - MED - Long - Mexico	Mexico
FBS308790	Venosan 5002 (Nylon) - Knee - 23-32Mmhg - Unisex - O/Toe - LGE - Long - Mexico	Mexico
FBS308800	Venosan 5002 (Nylon) - Knee - 23-32Mmhg - Unisex - O/Toe - XLGE - Long - Mexico	Mexico
FBS308810	Venosan 5002 (Nylon) - Knee - 23-32Mmhg - Unisex - C/Toe - SML - Long - Mexico	Mexico
FBS308820	Venosan 5002 (Nylon) - Knee - 23-32Mmhg - Unisex - C/Toe - MED - Long - Mexico	Mexico
FBS308830	Venosan 5002 (Nylon) - Knee - 23-32Mmhg - Unisex - C/Toe - LGE - Long - Mexico	Mexico
FBS308840	Venosan 5002 (Nylon) - Knee - 23-32Mmhg - Unisex - C/Toe - XLGE - Long - Mexico	Mexico
PTS645100	Venosan 5001 (Cotton) - Thigh W/ Belt - 18-22Mmhg - Unisex - O/Toe (Left) - SML - Long - Mexico	Mexico
PTS645110	Venosan 5001 (Cotton) - Thigh W/ Belt - 18-22Mmhg - Unisex - O/Toe (Left) - MED - Long - Mexico	Mexico
PTS645120	Venosan 5001 (Cotton) - Thigh W/ Belt - 18-22Mmhg - Unisex - O/Toe (Left) - LGE - Long - Mexico	Mexico
PTS645130	Venosan 5001 (Cotton) - Thigh W/ Belt - 18-22Mmhg - Unisex - O/Toe (Left) - XLGE - Long - Mexico	Mexico
PTS645140	Venosan 5001 (Cotton) - Thigh W/ Belt - 18-22Mmhg - Unisex - O/Toe (Right) - SML - Long - Mexico	Mexico
PTS645150	Venosan 5001 (Cotton) - Thigh W/ Belt - 18-22Mmhg - Unisex - O/Toe (Right) - MED - Long - Mexico	Mexico
PTS645160	Venosan 5001 (Cotton) - Thigh W/ Belt - 18-22Mmhg - Unisex - O/Toe (Right) - LGE - Long - Mexico	Mexico
PTS645170	Venosan 5001 (Cotton) - Thigh W/ Belt - 18-22Mmhg - Unisex - O/Toe (Right) - XLGE - Long - Mexico	Mexico

Code	Display Name	Colour
PTS645180	Venosan 5002 (Cotton) - Thigh W/ Belt - 23-32Mmhg - Unisex - O/Toe (Left) - SML - Long - Mexico	Mexico
PTS645190	Venosan 5002 (Cotton) - Thigh W/ Belt - 23-32Mmhg - Unisex - O/Toe (Left) - MED - Long - Mexico	Mexico
PTS645200	Venosan 5002 (Cotton) - Thigh W/ Belt - 23-32Mmhg - Unisex - O/Toe (Left) - LGE - Long - Mexico	Mexico
PTS645210	Venosan 5002 (Cotton) - Thigh W/ Belt - 23-32Mmhg - Unisex - O/Toe (Left) - XLGE - Long - Mexico	Mexico
PTS645220	Venosan 5002 (Cotton) - Thigh W/ Belt - 23-32Mmhg - Unisex - O/Toe (Right) - SML - Long - Mexico	Mexico
PTS645230	Venosan 5002 (Cotton) - Thigh W/ Belt - 23-32Mmhg - Unisex - O/Toe (Right) - MED - Long - Mexico	Mexico
PTS645240	Venosan 5002 (Cotton) - Thigh W/ Belt - 23-32Mmhg - Unisex - O/Toe (Right) - LGE - Long - Mexico	Mexico
PTS645250	Venosan 5002 (Cotton) - Thigh W/ Belt - 23-32Mmhg - Unisex - O/Toe (Right) - XLGE - Long - Mexico	Mexico