

Aspire Montana X Maxi Lift Recline Chair



Created for comfort and bariatric support, the **Montana X Maxi** features a zip-access adjustable backrest for personalised postural adjustments. Its sleek design fits effortlessly into modern homes, blending clinical function with elegant style.

✔ Contract approved

- **Trusted Bariatric Design:** 220kg high safe working load and generous seating surface
- **Refined, Contemporary Design:** Larger, plush armrests provide extra comfort while maintaining a streamlined look that suits home environments
- **Integrated Footrest Extension:** Reduce fatigue and lower limb swelling by improving healthy circulation
- **Protective Covers for High-wear Areas:** Removable headrest and armrest covers keep your chair clean, stylish, and long-lasting
- **Commercial-grade Seat Cushioning:** Enhances comfort and durability, providing lasting support
- **Promotes Independent Daily Living:** Simple 2-button handset assists users with cognitive or visual challenges

Warranty Information

- 2 Year Limited Warranty

Code	Display Name	Total Height	Safe Working Load	Material	Colour	Seat Depth	Seat Width	Backrest Height
CHP227401	Aspire Montana X Maxi - Golden Grain	1060mm	220kg	Microsuede	Golden Grain	530mm	640mm	670mm