

Avantgarde (CV, CS, CLT)



- Elegant design with consistent lightweight construction
- With its diversity of adjustment possibilities the Avantgarde is a great all-rounder
- Compact design provides maneuverability
- Fully customisable with a wide range of accessories available



ndis

Item Code:

05_122203108_0105_1_2

Wheelchair Manual –

Composite:05_417_0105_1_2

Click [here](#) for more information

ottobock.

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Safe Working Load	Unit Weight
MWC427115	Avantgarde Teen 2 (as per script)	240 - 400 mm	220 - 360 mm	175 - 400 mm	90 kg	10.7 kg