

Avantgarde (CV, CS, CLT)



ottobock.

- Elegant design with consistent lightweight construction
- With its diversity of adjustment possibilities the Avantgarde is a great all-rounder
- Compact design provides maneuverability
- Fully customisable with a wide range of accessories available

 Contract approved

Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_122203108_0105_1_2

Wheelchair Manual –

Composite:05_417_0105_1_2

Click [here](#) for more information

Code	Display Name	Safe Working Load	Unit Weight	Seat Depth	Seat Width
MWC427115	Avantgarde Teen 2 (as per script)	90 kg	10.7 kg	240 - 400 mm	220 - 360 mm