

Aidacare Aspire A205 Aluminium Patient Lifter



The Aidacare Aspire A205 is a compact and versatile aluminium mobile lifter with powered legs that is ideal for home, community and facility environments.

- **Transportation or Storage:** Is quick and convenient with tool free disassembly
- **Ergonomic Handset:** The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- **Adjustable Height:** Three adjustable height settings provides the caregiver with the flexibility to adjust the height range to suit the specific lift and patient scenario
- **Powered Leg Function:** Opens and closes the legs for lifting stability or accessibility
- **Quick Release Yoke Attachment:** One touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety
- **Build Quality:** Strong powder coated aluminium frame with industry leading Linak electronics and warranty provision
- **Safety:** Critical features include emergency stop button, emergency lowering and antientrapment technology, with both patient and caregiver safety at heart
- **Intelligent Diagnostics:** State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- **Tente Casters:** 2 x 75 mm front castors for maximum under bed clearance and accessibility. 2 x 100 mm rear locking castors for safety and manoeuvrability
- Item Weight: 38.1 kg
- Tested and compliant to AS/NZS ISO 10535:2011 (the Australian standard for hoists)

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123604167_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Accessories

Click on the cart icon to display the range of optional accessories

Warranty Information

- Frame: 5 Years
- Electronics: 3 Years
- Battery and Handset: 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Lift Range	Leg Gait Range	Electric Leg
LSS390610	Aidacare Aspire A205 Patient Lifter	1238 mm	651 - 1173 mm	1040 mm	205 kg	Aluminium	38.1 kg	415 - 1800 mm	547 - 1015 mm	True