

Aspire ComfiMotion Freestanding Allure Headboard



Aspire ComfiMotion headboards and footboards blend style with function, helping to create a dignified and practical sleeping environment. Whether freestanding or fixed, each variant addresses specific clinical and aesthetic needs - improving comfort for users and making life easier for carers.

Freestanding headboards stand independently behind the bed and can be moved or removed as needed. Contemporary in design, the headboard remains stationary while the bed is raised or lowered, not impacting any wall decoration above the headboard.

- Not attached to the bed - offers full positioning flexibility
- Quick to install or remove without tools
- Softened surface that supports pillows and protects the head
- Enhances room style without limiting bed function
- Supports high-needs equipment integration in care environments
- Available in two modern styles, both at 1140mm high and 110mm deep

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Compliance

IEC 60601-2-52:2009

Code	Display Name
BEA830120	Aspire ComfiMotion Freestanding Headboard - Allure - Portobello Grey - Long Single
BEA830130	Aspire ComfiMotion Freestanding Headboard - Allure - Portobello Grey - King Single
BEA830132	Aspire ComfiMotion Freestanding Headboard - Allure - Portobello Grey - Long Double
BEA830134	Aspire ComfiMotion Freestanding Headboard - Allure - Portobello Grey - Queen
BEA830136	Aspire ComfiMotion Freestanding Headboard - Allure - Portobello Grey - King

