

Softform Maxi Glide Mattress



- 2-surface “Glide” system reduces skin friction and shear whilst allowing the mattress to conform to changes in the bed functions
- Foam castellations on the upper surface of the mattress
- Two-way stretch PU waterproof cover
- Welded seams with full length flap conceals the zip
- Toughened PU coated base

Warranty Information

- 4 Years



 Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_043306003_0103_1_2

Power Bed –

Composite:05_423_0103_1_2

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
PTM636380	Softform Maxi Glide Mattress	1970 mm	880 mm	152 mm	248 kg