

Jobst Active Wear Compression Socks



JOBST ActiveWear offers cushioned compression socks designed for superior moisture management and exceptional softness. Featuring a wide comfort knee band, they are easy to put on and stay securely in place all day long. With a moderate 15-20 mmHg graduated compression, they are perfect for maintaining better circulation, making them the ideal choice for any active lifestyle.

- Dri-release yarns ensure dry and comfortable feet, resembling cotton.
- Wide comfort knee band keeps the sock in place.
- Cushioned foot with seamless toe provides all-day comfort.
- Recommended for tired, achy legs, swelling, discomfort from prolonged standing or sitting, preventing varicose veins and swelling during pregnancy, medical issues related to vein health, reducing swelling, and preventing DVT during travel.

Fabric Content: 60% Polyester, 20% Nylon, 11% Cotton, 9% Spandex



 Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name	Colour
PTS642700	Jobst ActiveWear - Knee High - 15-20mmHg - Closed Toe - White - Small	White
PTS642701	Jobst ActiveWear - Knee High - 15-20mmHg - Closed Toe - Black - Small	Black
PTS642702	Jobst ActiveWear - Knee High - 15-20mmHg - Closed Toe - White - Medium	White
PTS642703	Jobst ActiveWear - Knee High - 15-20mmHg - Closed Toe - Black - Medium	Black
PTS642704	Jobst ActiveWear - Knee High - 15-20mmHg - Closed Toe - White - Large	White
PTS642705	Jobst ActiveWear - Knee High - 15-20mmHg - Closed Toe - Black - Large	Black

Code	Display Name	Colour
PTS642706	Jobst ActiveWear - Knee High - 15-20mmHg - Closed Toe - White - XLarge	White
PTS642707	Jobst ActiveWear - Knee High - 15-20mmHg - Closed Toe - Black - XLarge	Black