

## Oppo 1033 knee support patella stabilizer



- Neoprene circle pad around the patella provides protection and stabilisation to prevent patella displacement
- Four neoprene straps provide support directly to patella area
- Two medial and lateral spiral stays provide additional support and stability for the knee
- Helps relieve pain from tendinitis and arthritis
- Retains body heat to warm knee joint and increases blood circulation
- Slips on with straps and Stays
- Note: Measure around the mid-section of the knee, while in a standing position.

✔ Contract approved

### Support @ Home

**Function:** Lower limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Knee Circumference
FBS301521	Oppo 1033 knee support patella stabilizer - Sml	33.7 - 37.5 cm
FBS301521M	Oppo 1033 knee support patella stabilizer - Med	37.5 - 40.6 cm
FBS301521L	Oppo 1033 knee support patella stabilizer - Lge	40.6 - 45.1 cm
FBS301521XL	Oppo 1033 knee support patella stabilizer - XL	45.1 - 50.2 cm