

Aspire Autoclavable Wash Bowls



Our range of Wash Bowls are fully autoclavable and made from premium polypropylene. They will last for many years in any healthcare environment subject to the rigours of cleaning, sanitation and infection control requirements.

- Heavy duty and reusable
- Designed for use in healthcare environments such as surgeries, acute care, general ward, medical centres, and aged care facilities
- Autoclavable polypropylene that can be cleaned using autoclaves and sanitisers
- Pastel green colour for easy identification
- *Length measurement refers to the diameter



Compliance

No applicable compliance standard

Warranty Information

- 1 Year



| Code | Display Name | RRP | Box Quantity | Total Length | Capacity |
|-----------|-------------------------|------|--------------|--------------|----------|
| WAT769800 | Wash Bowl - Small | 10.8 | - | 140 mm* | 600 ml |
| WAT770000 | Wash Bowl - Medium | 16.6 | - | 210 mm* | 1500 ml |
| WAT761200 | Wash Bowl - Large | 28.3 | 1 | 305 mm* | 2400 ml |
| WAT769600 | Wash Bowl - Extra Large | 37.4 | 1 | 345 mm* | 5500 ml |