

NuboMed A2310-S Medication Auxiliary Cabinet



NuboMed

- **Suitability:** Dispersed satellite inpatient area medication management and storage
- **User Interaction:** 1080p ten-point touch capacitive screen
- **Cabinet Material:** High-strength aluminium alloy and metal with anodised surface
- **Intelligent Lock Control:** Authorised user set up for cabinet lock access, allows individual set up for each bin
- **Cabinet Door:** Transparent tempered safety glass door
- **Drawer Layers:** 7 layers with 4 medicine bins per layer
- **Speaker:** Voice prompt with automatic error warning
- **Intelligent Monitoring:** Intelligent real-time quantity and specification monitoring, accommodates oral, ampoule, vial, boxed and blister pack medicines
- **Intelligent Display:** Each bin is equipped with a 3.5" colour LED display for drug identification, name, specifications, quantity and operating instructions
- **Anti-theft Detection:** Featuring an automatic alarm
- **Emergency Access:** Dual mechanical locks ensure that the medicine can still be accessed without power
- **Response Time:** <0.5 seconds detection
- **Ideal For:** Clinical Ward and Central Pharmacy



Warranty Information

- 2 Years