

3D Sideboard



- Refined sideboard that blends functional design with aesthetic appeal
- Suited to a hall way, living or private area
- Featuring three spacious drawers and cupboards to store magazines, games and personal items
- Slimline curved handles add a distinct mid-century feel
- Generous surface area to display decorative items and treasured mementos
- Dimensions: 900mm (H) x 1500mm (L) x 450mm (W)
- Pictured in Milano Walnut, colour/material to be advised at time of order



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
FUL053120	3D Sideboard	1500mm	450mm	900mm