

AC3 King Single Bed 3-Section - 1090mm w (220kg)



HIRE ME

The Aidacare AC3 bed is a three function bed tailored to the Aged Care environment. Designed in Australia to meet International Standards, the premium frame and Linak electronics will ensure long term reliability and safety. Multiple functions and accessories combine to exceed the comfort, support and safety expectations for a diverse range of users.



- Leading Quality Frame & Electronics - Australian designed durable steel frame with European Linak electronics and handset allows and handset allows for market leading warranty provision
- Integrated Bed Extension - Provides an additional 195mm length to cater for all users
- Central Locking - 125mm twin castors with full width foot pedal
- Wall Protection System - A wall locator bar A wall locator bar facilitates bed positioning facilitates bed positioning and helps prevent wall and helps prevent wall damage
- Cable Management - Ensures cables are Ensures cables are kept safely out of the kept safely out of the way for residents and way for residents and care staff. Easy access service compartment for service compartment for electronics
- Removable Laminate Bed Ends - In a contemporary beech finish. Any finish available to order

*** Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

Code	Display Name
H_BEB045905	AC3 King Single Bed 3-Section - 1090mm w (220kg)