

Two Handled Mug



- Designed to make drinking safer and easier for people with hand tremors
- Supplied with two lids, one for anti-spill and one with a drinking spout and two handles to assist with control
- Sold as a pair (pack includes 2 cups, 4 lids)
- Weight: 100 g

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Eating and drinking

Category: Low risk

Click [here](#) for more information

NDIS

Item Code:

05_150900319_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Capacity
DLK276520	Two Handled Mug	200 ml