

Kerry Pivot Frames



- Carer controlled frame to pivot and manoeuvre a patient during transfers
- Easily adjustable angle from 10 to 118 degrees off verticle
- Padded head section surrounding frame for user safety and comfort
- Allows a more comfortable user seating position and more control for the carer
- Safe Working Load: 225 kg

Warranty Information

- 2 years



Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load
LSS383600	Pivot Frame - Keyhole Sling Points	225 kg