

DermaSaver Knee Tube



- Soft and loose with ample airflow reaching the skin allowing for greater healing ability whilst offering protection
- Cushioned layers of MicroSpring Textile protect fragile skin from friction and pressure that can cause skin tears, bruising and abrasions
- Sizing: simply measure the circumference around the widest part of the calf. For borderline measurements, choose the larger size to ensure comfort

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Calf/Arm Circumference
PTA533320	DermaSaver - Knee Tube - X Sml	18-23cm
PTA533320S	DermaSaver - Knee Tube - Sml	23-30cm
PTA533320M	DermaSaver - Knee Tube - Med	30-38cm
PTA533320L	DermaSaver - Knee Tube - Lge	38-48cm
PTA533320XL	DermaSaver - Knee Tube - XL	48-58cm
PTA533321M	DermaSaver - Double Knee Tube - Med	30-38cm
PTA533321L	DermaSaver - Double Knee Tube - Lge	38-48cm

Code	Display Name	Calf/Arm Circumference
-------------	---------------------	-------------------------------