

Aspire Homecare Shower Chair



- Provides a complete seating support for users with reduced mobility during showering
- Durable high grip feet for stability and safety
- Adjustable legs to suit a wide range of users
- Raised armrests assist the user in transferring on and off
- Width between the arms: 465 mm

- ✔ Contracted DVA product
- ✔ Contract approved

Warranty Information

- 1 Year

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information



Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
BTS118025	Aspire Homecare Shower Chair	525 mm	750 - 880 mm	100 kg	340 mm	410 mm	440 - 570 mm	440 - 570 mm