

## Freya Premium Bedside Cabinet



- Bedside cabinet offering a contemporary take on traditional styling
- Sleek tapered legs to complement a variety of interior styles
- Features a lockable top drawer to store valuables and medication
- Generous cupboard storage for books, magazines and larger essentials
- Ample space on top for a lamp and smart devices
- Choose from an extensive range laminates or veneer stain options
- Pictured in Milano Walnut



### Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
BEA828450	Contemporary Bedside Cabinet	500mm	480mm	700mm
BEA826280	Freya Premium Bedside Cabinet - Timber	450 mm	500 mm	700 mm