

DONJOY 82915/6 RESPIFORM WRIST SUPPORT



- Dedicated for either right or left hand use
- Ergonomic, flexible lightweight design
- Anatomically designed, provides tailored individual fit
- Easy and simple application via adjustable straps and 'Stop Strap' system
- Incorporates two semi-rigid dorsal stays and an aluminum palmar stay for overall wrist immobilisation
- Advanced Ozone™ material, promotes air flow
- Breathable in general use, providing good skin care and around-the-clock comfort in different heat and humidity conditions Application:
- Wrist sprains and trauma
- Tendinopathy
- Pre and post-operative immobilization

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS305850	DONJOY 82915/6 RESPIFORM WRIST SUPPORT - LEFT - XS
FBS305860	DONJOY 82915/6 RESPIFORM WRIST SUPPORT - LEFT - SML
FBS305930	DONJOY 82915/6 RESPIFORM WRIST SUPPORT - RIGHT - LGE
FBS305940	DONJOY 82915/6 RESPIFORM WRIST SUPPORT - RIGHT - XL
FBS305870	DONJOY 82915/6 RESPIFORM WRIST SUPPORT - LEFT - MED
FBS305880	DONJOY 82915/6 RESPIFORM WRIST SUPPORT - LEFT - LGE
FBS305890	DONJOY 82915/6 RESPIFORM WRIST SUPPORT - LEFT - XL

Code	Display Name
FBS305900	DONJOY 82915/6 RESPIFORM WRIST SUPPORT - RIGHT - XS
FBS305910	DONJOY 82915/6 RESPIFORM WRIST SUPPORT - RIGHT - SML
FBS305920	DONJOY 82915/6 RESPIFORM WRIST SUPPORT - RIGHT - MED