

K Care URETHANE Commode Seat Padded - Open



- High quality urethane padded seat with an open front for personal cleaning and toileting access
- Available in blue
- Supplied with plastic seat runners for easy attachment and removal during cleaning of the seat and commode
- Seat Height: 40 mm

Warranty Information

- 1 Year



Contract approved

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width
BTC075700	Open Front Commode Seat - 445 mm Wide	430 mm	445 mm