

Cura - Spare Parts and Accessories

Accessories and spare parts for the Cura alarms product range.



✓ Contract approved

Support @ Home

Function: Signalling, alarming and localisation

Category: Prescribed

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width
BEA007460	Pull String Monitor	-	-
BEA007462	Wristband Transmitter Strap	255 mm	19 mm
BEA007463	Cura - Wireless Waterproof Call Button	-	-
BEA007466	Cura Alarms - Personal - ActiveCare Pager	-	-
BEA007550	Cura - ActiveCare Magnetic Reed Switch	-	-
BEA007549	Cura - ActiveCare Wireless Handset (includes wall bracket)	-	-
BEA007548	Cura - Nurse call pendant with 2m Cord (for BEA007409, BEA007410, BEA007428)	-	-

Code	Display Name	Total Length	Total Width
BEA007546	Cura - Monitor Soft Cover For Monitors with rubber strap (for BEA007410, BEA007428)	-	-
BEA007545	Cura - Wall Mounting Bracket (for BEA007428, BEA007410, BEA007466)	-	-
BEA007543	Cura - Caregiver Reset Key (for BEA007410, BEA007428)	-	-
BEA007457	Cura - EzyStand for iBeam motion sensor	-	-
BEA007482	Cura - Bed Monitor Plug Pack 6V (replacement)	-	-