

GloriaMed (MicroMesh) Pantyhose



- Fashionable and sheer
- Suits all standard body shapes
- 4 Sizes
- Short or long fittings
- Closed-toe
- Beige or black

Gloria Med general overview:

- Made in Italy, Gloria Med is one of the oldest and most respected manufacturers of compression legwear
- Conforms to all recognised international standards

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name	Colour
PTS641877BL0	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C /Toe - Sml -Long- Blk	Black
PTS641877BL1	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C /Toe - Med -Long- Blk	Black
PTS641877BL2	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C /Toe - Lge -Long- Blk	Black
PTS641877BL3	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C /Toe - XLge -Long- Blk	Black
PTS643940	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C/Toe - Sml -Short- Bge	Beige
PTS643950	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C/Toe - Sml -Short- Blk	Black
PTS643960	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C/Toe - Sml -Long- Bge	Beige

Code	Display Name	Colour
PTS643970	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C/Toe - Med -Short- Bge	Beige
PTS643980	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C/Toe - Med -Short- Blk	Black
PTS643990	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C/Toe - Med -Long- Bge	Beige
PTS644000	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C/Toe - Lge -Short- Bge	Beige
PTS644010	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C/Toe - Lge -Short- Blk	Black
PTS644020	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C/Toe - Lge -Long- Bge	Beige
PTS644030	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C/Toe - XLge -Short- Bge	Beige
PTS644040	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C/Toe - XLge -Short- Blk	Black
PTS644050	GloriaMed - MicroMesh P/hose - 18-21mmHg - Female - C/Toe - XLge -Long- Bge	Beige
PTS644540	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - Sml -Short- Bge	Beige
PTS644550	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - Sml -Short- Blk	Black
PTS644560	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - Sml -Long- Bge	Beige
PTS644570	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - Sml -Long- Blk	Black
PTS644580	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - Med -Short- Bge	Beige
PTS644590	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - Med -Short- Blk	Black
PTS644600	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - Med -Long- Bge	Beige
PTS644610	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - Med -Long- Blk	Black

Code	Display Name	Colour
PTS644620	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - Lge -Short- Bge	Beige
PTS644630	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - Lge -Short- Blk	Black
PTS644640	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - Lge -Long- Bge	Beige
PTS644650	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - Lge -Long- Blk	Black
PTS644660	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - XLge -Short- Bge	Beige
PTS644670	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - XLge -Short- Blk	Black
PTS644680	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - XLge -Long- Bge	Beige
PTS644690	GloriaMed - MicroMesh P/hose - 23-32mmHg - Female - C/Toe - XLge -Long- Blk	Black