

Aspire Elevating Legrest



- Assists in leg and foot positioning or elevation, most commonly to treat oedema
- Legrest is angle adjustable to six positions
- Footplate is height and angle adjustable
- Fits Aspire Shower Commodes

Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS 3973:2009



Code	Display Name
BTC067110S	Elevating Legrest - Tilt in Space - Left
BTC067115S	Elevating Legrest - Tilt in Space - Right
BTC067108M	Elevating Legrest - Non-Tilt Commodes - LEFT
BTC067109M	Elevating Legrest - Non-Tilt Commodes - RIGHT
BTC067110M	Elevating Legrest - Tilt in Space - LEFT
BTC067115M	Elevating Legrest - Tilt in Space - RIGHT