

## Jay - ZIP Backrest



**JAY®**

- Designed specifically to provide posterior and mild lateral thoracic support for active kids aged 2 to 15 years of age
- The contour shapes vary by size, designed to accommodate the unique shapes of clients of each width. The Zip's extra soft foam provides superior immersion for children's lower weights
- Versatile hardware offers 3" of depth adjustment and 1" of width adjustment to accommodate kids' rapid growth. Its 40° of angle adjustability helps to improve visual orientation and function
- The Zip features an innovative, dual-cover system. The inner cover is water-resistant and easy to wipe off and clean, while its Aqua-guard zipper and anti-wicking thread protect the foam from fluids. The machine-washable, stretch outer cover utilizes naturally anti-microbial, silver-impregnated X-static® fiber
- All models include: X-static® silver thread stretch cover and universal mounting hardware

### Warranty Information

- 2 Years Backrest
- 6 Months Cover

 Contract approved

### Support @ Home

**Function:** Accessories for wheelchair seating

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_181003192\_0105\_1\_2

**Wheelchair Manual – Composite:**

**Wheelchair Powered –**

**Composite:**05\_417\_0105\_1\_2

05\_419\_0105\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Safe Working Load	Seat Width	Backrest Height	Backrest Width
PTB533844	JAY ZIP - 10 x 10in Tall	75 kg	9 - 11 in	10 in	10 in

Code	Display Name	Safe Working Load	Seat Width	Backrest Height	Backrest Width
PTC585194	JAY ZIP - 12 x 8in Tall	75 kg	11 - 13 in	8 in	12 in
PTB533847	JAY ZIP - 12 x 10in	75 kg	11 - 13 in	10 in	12 in
PTB533845	JAY ZIP - 12 x 12in Tall	75 kg	11 - 13 in	12 in	12 in
PTB533848	JAY ZIP - 14 x 12in X Tall	75 kg	13 - 15 in	12 in	14 in