

## Kerry Full Body Sling - Polyester



- Attaches to hoist to enable a carer to lift and manoeuvre a less mobile user for toileting
- Provides excellent back and head support for patients requiring support during transfers
- Divided leg is ideal for toileting patients who require full body support
- Quick and easy to fit
- Colour coded attachment and positioning straps
- Machine washable
- FIM Scores: 1 and 2



**ndis**

**Item Code:**

05\_123603165\_0105\_1\_2

**Sling – Composite:**

05\_421\_0105\_1\_2

Click [here](#) for more information

### Warranty Information

- 1 Year

**Compliance**

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS382400	Full Body Sling - Medium	250 kg	Polyester
LSS382500	Full Body Sling - Large	250 kg	Polyester
LSS382550	Full Body Sling - Extra Large	250 kg	Polyester